

# Medela breast pump component compatibility chart

To ensure proper performance and avoid breast pump damage, use only the components designed to work with your breast pump.

Breast Pumps		Compatible Connector	Compatible Tubing	Part Number
	<b>Pump In Style®</b> WITH <b>maxFLOW™</b>	 OR 	 Single Y- tubing : All three adaptor ends are narrow & round.	<b>Connector</b> 101041267 101038233 <b>Tubing</b> 101040485
	<b>Freestyle Flex™</b>	 OR 	 Single Y- tubing : two narrow round ends and 1 larger triangle end.	<b>Connector</b> 101038223 101041267 <b>Tubing</b> 101038234
	<b>Sonata®</b>		 Single Y- tubing : two narrow round ends and 1 larger round end.	<b>Connector</b> 68054 <b>Tubing</b> 101038540
	<b>Harmony®</b>	 OR 	NO TUBING	<b>Connector</b> 101034753 101034864
	<b>Swing®</b>		 One single 30" tube: with one round end and one triangle end.	<b>Connector</b> 101034753 <b>Tubing</b> 101034768
	<b>Symphony®</b>		 Two separate 42" tubes: with one round end and one triangle end.	<b>Connector</b> 87071 <b>Tubing</b> 101034773
	<b>Pump In Style® Advanced</b>	 OR 	 Two separate tubes: with one round end and one soft tubing end.	<b>Connector</b> 101039190 87071 <b>Tubing</b> 101033078
	<b>Freestyle™</b>		 Single Y- tubing : All three adaptor ends are thick & round.	<b>Connector</b> 67061 <b>Tubing</b> 101033077